



# Medicinal Chemistry

## The first stage

### College of Dentistry



By

Assistant Lecturer

***HAZIM SAAD JABBAR AL-MALIKI***

University of Basra - College of Dentistry - Department  
of Basic Sciences

MSc. Medicinal Biochemistry from Malaysian National  
University (UKM) – Malaysia

Email:- [hazim.saad@uobasrah.edu.iq](mailto:hazim.saad@uobasrah.edu.iq) and [hazzmrn@yahoo.com](mailto:hazzmrn@yahoo.com)



**QUANTITATIVE  
VOLUMETRIC  
ANALYSIS**

# **CONTENTS:-**

**1**

**INTRODUCTION**

**2**

**What is the meaning of Titration?**

**3**

**A reagent , A primary standard solution and Types of titrations**

**4**

**Acid-base, Precipitation, Redox and Complexometric titration**

**5**

**Calibration of Hydrochloric Acid**




# INTRODUCTION

## Volumetric analysis

**Volumetric analysis is a general term for a method in quantitative chemical analysis in which the amount of a substance is determined by the measurement of the volume that the substance occupies. It is commonly used to determine the unknown concentration of a known reactant. Volumetric analysis is often referred to as titration.**

## What is the meaning of Titration?

**Titration is a common laboratory method of quantitative chemical analysis that is used to determine the unknown concentration of a known reactant. Because volume measurements play a key role in titration, it is also known as volumetric analysis.**




A reagent, called the titrant or titrator,[1] of a known concentration (a standard solution) and volume is used to react with a solution of the analyte or titrant,[2] whose concentration is not known. Using a calibrated burette or chemistry pipetting syringe to add the titrant.

A **primary standard solution** is a highly purified compound that serve as a reference material in all volumetric titrimetric methods.

## Types of titrations

- 1- Acid-base titration
- 2- Precipitation titration
- 3- Redox titration
- 4- Complexometric titration

# Acid-base titration



# Calculations of volumetric analysis

**Standard solution** is one, which contains a known weight of the reagent in a definite volume of the solution.

**Molar solution** is one, which contains 1 gm molecular weight of the reagent per liter of solution.

$$M = \frac{\text{Weight}}{\text{M.Wt}} * \frac{1000}{\text{Volume(mL)}}$$

**Normal solution** is one that contains 1gm equivalent weight per liter of solution.

$$N = \frac{\text{Weight}}{\text{eq.wt}} * \frac{1000}{\text{Volume(mL)}}$$



## Equivalent weights

(1) Equivalent weight in neutralization reactions.

**The equivalent weight of acid is that weight of it which contains one-gram atom of replaceable hydrogen.**

**Ex: equivalent weight of  $\text{H}_2\text{SO}_4 = \text{M.Wt } \text{H}_2\text{SO}_4 / 2$**

**equivalent weight of  $\text{H}_3\text{PO}_4 = \text{M.Wt } \text{H}_3\text{PO}_4 / 3$**

$$\text{eq.wt acid} = \frac{\text{M.Wt acid}}{\text{No. of active H}^+}$$

**The equivalent weight of Base is that weight of it which contains one replaceable hydroxyl group.**

**Ex: equivalent weight of  $\text{NaOH} = \text{M.Wt } \text{NaOH} / 1$**



# Experiment (5)

## Calibration of Hydrochloric Acid

### A. Preparation hydrochloric acid solution (0.1N):-

1- Calculate the normality of the concentrated HCl :

$$N1 = \frac{\text{SP.} * \text{Wt.}\% * 1000}{\text{Eq. Wt.}} \quad \text{----- (1)}$$

Where :-

N1 = hydrochloric acid calibration Centre

Sp. = Specific weight of acid (acid density)

Wt. % = The percentage of the weight of hydrochloric acid Center

$$N = \frac{1.19 * (37/100) * 1000}{36.5} = 12.0630$$



## 2- To prepare (500mL) of 0.1N HCl

Preparation titrated hydrochloric acid ( $N_2 = 0.1$ ) and size ( $V_2 = 500$  ML) Use the following formula to calculate the size of the acid Centre ( $V_1$ ), which we take to prepare the acid diluted:

$$N_1 * V_1 = N_2 * V_2 \quad \text{-----}(2)$$

conc. HCl    dil. HCl

$$V_1 = \frac{N_2 * V_2}{N_1} \quad \text{-----}(3)$$


$$12.0630 \times V_1 = 500 \times 0.1$$
$$V_1 = 4.1449 \text{ mL}$$

## B: Preparation of (0.1N) sodium carbonate (Na<sub>2</sub>CO<sub>3</sub>):

To calculate the weight of sodium carbonate needed to prepare a titration solution (0.1), the following equation is used:

$$N = \frac{\text{Weight}}{\text{equivalent weight}} * \frac{1000}{\text{Volume(mL)}}$$

$$\text{Wt} = \frac{N * V * \text{Eq.Wt}}{1000}$$



Where:

WT = weight of sodium carbonate dissolved in a volume of water

N = Standard of Sodium Carbonate to be prepared = 0.1

V = volume of the solution to be prepared (in ml)

Eq.Wt = equivalent weight to sodium carbonate



## Equipment:-

1) burette of 50 ml. 2) 10 ml pipette. 3) 500ML Beaker. 4) 500 ML volumetric vial, 250 ML volumetric flask and conical flask. 5) funnel. 6) Spatule. 7) cylinder.

## Procedure: -

### A. Standardization of HCl solution with standard solution of $\text{Na}_2\text{CO}_3$

1-Clean the burette and rinse with HCl solution.

2-Fill the burette with HCl.

3-Pipet 10 ml of standard solution ( $\text{Na}_2\text{CO}_3$ ) into a 250mL conical flask. Add 3 drops of **Methyl orange indicator**.

4-Titrate by adding HCl drop wise until the solution just begins to change from yellow to red.


5-Repeat the titration a few times until you get approximate results. Take the average of the results and subtract 0.05mL. (This result represents the volume of extra drop which change the color of indicator.

6-Calculate the normality of HCl :


$$N_{\text{acid}} \times V_{\text{acid}} = N_{\text{base}} \times V_{\text{base}}$$

7- Make label on your bottle containing your name, date of preparation and concentration of acid after standardization.

8- The equation of reaction



## A Presentation On Precipitation Titration






## B. Analysis of sodium carbonate $\text{Na}_2\text{CO}_3$

- 1-Clean the burette and rinse with standardized HCl solution and then fill it with the acid.
- 2-Pipet 10 ml unknown solution ( $\text{Na}_2\text{CO}_3$ ) into a 250mL conical flask. Add 2 drops of **phenolphthalein indicator** the solution will be pink.
- 3-Titrate by adding HCl drop wise until the solution just begins to change its color from pink to **colorless** this data will be ( $V_1$ ).
- 4- Add 1-2 drops of **Methyl orange indicator** to the above solution which became yellow then complete the titration until the color of the solution became pale **orange (onion)**, this data will be ( $V_2$ ).
- 5-Repeat the titration a twice time until you gets approximate results. Take the average of the results and subtract 0.05mL. (This result represents the volume of extra drop which change the color of indicator.
- 6- Make a table as bellow: -

Sq		First titration	Second titration	Third titration	Average of titrations
1	Titration with ph.ph	V1	V1	V1	V1 (av.)
2	Titration with M.O	V2	V2	V2	V2 (av.)







$$V1(av.) = \frac{V1+V1+V1}{3}$$



$$N_{acid} \times V_{(tot.) acid (from burette)} = N_{base} \times V_{base}$$

7- The equation of reaction :




# Methyl Orange Indicator

Acidic pH



Red - Orange



Neutral pH




Yellow



Alkaline pH



Yellow



# Phenolphthalein indicator

pH = 0-3




Orange



pH = 4-7



Colorless



pH = 8-10




Pink



pH = 11-14



Colorless





**Thank You For Listening**