



Medicinal Chemistry

The first stage

College of Dentistry



By

Assistant Lecturer

HAZIM SAAD JABBAR AL-MALIKI

University of Basra - College of Dentistry - Department
of Basic Sciences

MSc. Medicinal Biochemistry from Malaysian National
University (UKM) – Malaysia

Email:- hazim.saad@uobasrah.edu.iq and hazzmrn@yahoo.com

Preparation of different types of solutions

Experiment (6)

Standardization of sodium hydroxide

NaOH solution with standard solution of
hydrochloric acid HCl



Important requirements for primary standard compounds are:

- **Extremely pure**
- **Highly stable**
- **Anhydrous**
- **Can be weighed easily**
- **Should be ready to use and available**
- **Should be preferably non toxic**
- **Should not be expensive**

Objectives

- **Standardize a sodium hydroxide solution**
- **Determine the Molarity or Normality of an unknown hydrochloric acid solution**
- **Understand the use of indicators in titrations**
- **Learn proper pipetting technique**
- **Learn to titrate a strong acid with a strong base.**


Preparation of standard solution of Na_2CO_3 (0.1N):

- 1- Weigh out accurately 1.325gm of Na_2CO_3 .
- 2- Dissolve in small quantity of distilled water and transfer quantitatively to 250ml measuring flask.
- 3- Complete to the mark and shake well.
- 4- Calculate the exact normality of Na_2CO_3 solution.



$$N = \frac{\text{Weight}}{\text{eq.wt}} * \frac{1000}{\text{Volume(mL)}}$$

Weight required = Normality x eq.wt. x volume in (1000 ml (liter)).



(I) Determination of the normality of hydrochloric acid by a standard solution of sodium carbonate (0.1N).

Sodium carbonate reacts with hydrochloric acid according to the following equation:




When one equivalent of HCl is added to the carbonate it is transformed into bicarbonates.



pH= 11.5 (alkaline) to 8.3.

HCl is added to the solution of bicarbonate



pH = 8.3 to 3.8

Calculations:

- **In case of M.O.**

$$N_1 V_1(\text{HCl}) = N_2 V_2(\text{Na}_2\text{CO}_3)$$

Or


$$N_1 = N_2 V_2 / V_1$$

- **In case of ph.ph .**

$$N_1 2V_1 = N_2 V_2$$

Or

$$N_1 = N_2 V_2 / 2V_1$$



(II) Determination of the strength and normality of sodium hydroxide solution by a standard solution of hydrochloric acid

HCl reacts with sodium hydroxide according to the following equation:



Calculations:

In both cases of M.O. and ph.ph.

use the relation:

$$N = \frac{\text{Weight}}{\text{eq.wt}} * \frac{1000}{\text{Volume(mL)}}$$

$$N_1 V_1(\text{HCl}) = N_2 V_2(\text{NaOH})$$

In order to deduce the normality of NaOH



Thank You For Listening