

Preparation of different types of solutions

Experiment (7)

An analysis of a mixture of Na_2CO_3 and NaOH using two indicators and a standard HCl solution

Theory:

1- When a known volume of the mixture is titrated with HCl in presence of ph. ph., the acid reacts with all the sodium hydroxide and with only half of the carbonate.

$$V1 = \text{all hydroxide} + 1/2 \text{ the carbonate}$$

2- When a known volume of the mixture is titrated with HCl in presence of M.O., the acid reacts with all the hydroxide and all the carbonate.


$$V2 = \text{all hydroxide} + \text{all carbonate}$$

$$\text{Volume of HCl} = 1/2 \text{ carbonate} = V2 - V1 = V \text{ ml}$$

$$\text{Volume of HCl} = \text{all carbonate} = 2V \text{ ml}$$

$$\text{Volume of HCl} = \text{NaOH} = V2 - 2V \text{ ml}$$

General interaction



Materials and method:

1- Fill the burette with 0.1N (HCl) solution.

2- Transfer with a pipette 10 ml of the mixture (NaOH + Na₂CO₃) to a conical flask, and add one or two drops of ph. ph.

3- Add the acid from the burette till the solution becomes colorless.

4- Repeat the experiment twice or three times and tabulate your results.

- The volume of acid in the case is V_x and is equivalent to all the hydroxide and half the carbonate (NaOH + $\frac{1}{2}$ Na₂CO₃).

5- Repeat the experiment with methyl orange (M.O) until the color of the solution is changed to faint red.

6- Repeat the experiment twice or three times and tabulate your results.

- The volume of acid in the case is V_y and is equivalent to all the hydroxide and all the carbonate (NaOH + Na₂CO₃), which represents the volume of the acid, which is equivalent to ($\frac{1}{2}$ Na₂CO₃).

7- Calculate the strength of the sodium hydroxide and the sodium carbonate in the mixture.

Calculations:

- **In the case of Na_2CO_3 :**

- Calculate the concentration of (Na_2CO_3).

- Calculate the volume of acid (HCl) that is equivalent to (Na_2CO_3)

$$V_1 = 2V_y$$

$$N_{\text{HCl}} \times V_{\text{HCl}} = N_{\text{Na}_2\text{CO}_3} \times V_{\text{Na}_2\text{CO}_3}$$

- **In the case of NaOH:**

- Calculate the volume of acid equivalent (NaOH).

$$V_2 = V_x - V_y$$

$$V_2 = V_x (1/2 \text{Na}_2\text{CO}_3 + \text{NaOH}) - V_y (1/2 \text{Na}_2\text{CO}_3)$$

$$N_{\text{HCl}} \times V_{\text{HCl}} = N_{\text{NaOH}} \times V_{\text{NaOH}}$$

Experiment (8)

Precipitation titration


Determination of chloride ion by Mohr method

Precipitation titration: is titration depend upon the combination of ions to form a simple precipitate. Mohr method is a method depend upon formation a colored precipitate for the determination of chloride ion.

Chloride ion, is reacted with silver nitrate solution to form AgCl precipitate.



A small quantity of potassium chromate (K_2CrO_4) solution is added to serve as indicator. The first excess of titrant results in the formation of a red silver chromate precipitate which signal the end point.



Procedure:

- 1-Clean the burette and fill it with silver nitrate (0.1N).
- 2-Pipet 10mL of chloride ion solution into 250mL conical flask, add 5drops of potassium chromate.
- 3-Titrate chloride solution against silver nitrate until arrive to equivalent point(the point in which the number of moles of AgNO_3 equal to the number of moles of chloride ion)[Notice a white precipitate in the yellow solution].

After this point the excess of AgNO_3 will react with potassium chromate leading to formation of red precipitate Ag_2CrO_4 (the end point). The difference between equivalent and end point is the volume of AgNO_3 reacted with the indicator.

4- Repeat the titration and calculate the average and subtract the indicator blank (equal to 0.1mL). Calculate the normality.

$$N_{\text{Cl}^-} \times V_{\text{Cl}^-} = N_{\text{Ag}^+} \times V_{\text{Ag}^+}$$

$$\text{Concentration of } [\text{Cl}^-] (\text{ppm}) = N_{\text{Cl}^-} \times \text{eq.wt} \times 1000$$