



Medicinal Chemistry

The first stage

College of Dentistry



By

Assistant Lecturer

HAZIM SAAD JABBAR AL-MALIKI

University of Basra - College of Dentistry - Department
of Basic Sciences

MSc. Medicinal Biochemistry from Malaysian National
University (UKM) – Malaysia

Email:- hazim.saad@uobasrah.edu.iq and hazzmrn@yahoo.com

Preparation of different types of solutions

Experiment (7)

An analysis of a mixture of Na_2CO_3 and NaOH using two indicators and a standard HCl solution



Theory:

1- When a known volume of the mixture is titrated with HCl in presence of ph. ph., the acid reacts with all the sodium hydroxide and with only half of the carbonate.

$$\mathbf{V1 = all\ hydroxide + 1/2\ the\ carbonate}$$

2- When a known volume of the mixture is titrated with HCl in presence of M.O., the acid reacts with all the hydroxide and all the carbonate.


$$\mathbf{V2 = all\ hydroxide + all\ carbonate}$$

$$\mathbf{Volume\ of\ HCl = 1/2\ carbonate = V2 - V1 = V\ ml}$$

$$\mathbf{Volume\ of\ HCl = all\ carbonate = 2V\ ml}$$

$$\mathbf{Volume\ of\ HCl = NaOH = V2 - 2V\ ml}$$

General interaction




Calculations:

- In the case of Na_2CO_3 :
- Calculate the concentration of (Na_2CO_3).
- Calculate the volume of acid (HCl) that is equivalent to (Na_2CO_3)

$$V_1 = 2V_2$$

$$N_{\text{HCl}} \times V_{\text{HCl}} = N_{\text{Na}_2\text{CO}_3} \times V_{\text{Na}_2\text{CO}_3}$$



- **In the case of NaOH:**

- Calculate the volume of acid equivalent (NaOH).

$$V_2 = V_X - V_y$$

$$V_2 = V_X (1/2 \text{ Na}_2\text{CO}_3 + \text{NaOH}) - V_y (1/2 \text{ Na}_2\text{CO}_3)$$

$$N_{\text{HCl}} \times V_{\text{HCl}} = N_{\text{NaOH}} \times V_{\text{NaOH}}$$



Experiment (8)

Precipitation titration

Determination of chloride ion by Mohr method

Precipitation titration: is titration depend upon the combination of ions to form a simple precipitate. Mohr method is a method depend upon formation a colored precipitate for the determination of chloride ion.

Chloride ion, is reacted with silver nitrate solution to form AgCl precipitate.



A small quantity of potassium chromate (K_2CrO_4) solution is added to serve as indicator. The first excess of titrant results in the formation of a red silver chromate precipitate which signal the end point.



Calculate

$$N_{Cl^-} \times V_{Cl^-} = N_{Ag^+} \times V_{Ag^+}$$

$$\text{Concentration of } [Cl^-](\text{ppm}) = N_{Cl^-} \times \text{eq.wt} \times 1000$$



Thank You For Listening