

## **Laboratory instruments**

- 1- The balance of sensitive
- 2- Pipette
- 3- Conical flask
- 4- Burette
- 5- Volumetric flask
- 6- Beaker
- 7- Cylinder
- 8- Funnel
- 9- Watch class
- 10- Filter paper
- 11- Stand
- 12- Clamp

\* Types of Titration:

1- Acid-base Titration

2- Redox Titrations


3- Precipitation Titration

4- Complexometric Titration

### \*Experiment 1

Name of experiment: Titration 0.1N hydrochloric acid (HCl) with sodium hydroxide (NaOH)

Purpose of experiment: knew the concentration right for HCl.



In the titration of hydrochloric acid (HCl) with a base such as Sodium hydroxide (NaOH), 1 molecule of HCl will react with 1 molecule of NaOH to produce 1 molecule of the salt, Sodium Chloride (NaCl) and 1 molecule of water. The chemical equation allows us to calculate the concentration of a solution of HCl by titration with the base NaOH (where the concentration of NaOH is accurately known). Recall that at the equivalence point the amount of base added is chemically equal to the amount of acid present in the solution.

$$N_{\text{HCl}} V_{\text{HCl}} = N_{\text{NaOH}} V_{\text{NaOH}}$$

$\text{moles HCl} = \text{moles NaOH}$
--

\*Equivalent point: is a point in which the added titration is chemically equivalent to the analyte in the sample.

\***End point:** is a point where the indicator changes its colour.