

Organic Pharmaceutical Chemistry

^{3rd}Stage ^{1st} Semester

Lab No: 3



Assay of citric acid

Prepared by:

Assist. lecturar: Mustafa Mohammed Hussein Al-Hakiem

- B.Sc. Pharmacy 2013.

- M.Sc. Pharmaceutical Chemistry 2020.

-E. Mail:mustafa.hashim@uobasrah.edu.iq

Introduction:

Citric acid (C₆H₈O₇, m.wt.=192.123) is a white crystalline powder (or colorless crystals).

- It is very soluble in water and freely soluble in alcohol.
- It is a tri-basic acid, so its solutions are strongly acidic.
- It is available as the anhydrous form or monohydrate form. Assay of citric acid is based on the anhydrous form.





Structure of citric acid

Chemical principle:

Since citric acid has strong acid properties, it is titrated against a standard basic solution like 1N NaOH solution in an acid- base reaction.

 $H_2 - C - COOH + 3NaOH = H_2 - C - COONa + 3H_2O + H_2 - C - COONa + 3H_2O + H_2 - C - COOH + 3NaOH = H_2 - C - COONa + 3H_2O + H_2 - C - COONA + H_2 - C - C - COONA + H_2 - C - C - COONA + H_2 -$

Reaction between citric acid and sodium hydroxide solution

Uses of citric acid:

Uses in medicine and pharmacy:

Citric acid is used widely in effervescent salts, for their effervescence on CO₂ produced from the reaction between citric acid and sodium bicarbonate

Citric acid is also used to dissolve renal stones.

Sodium citrate used in blood transfusion and bacteriology for the prevention of blood clotting.



EFFERVESCENT POWDER

Relief from painful burning symptoms of urinary tract infections

URINARY ALKALINIZER

Each 4 g sachet contains: sodium bicarbonate 1.76 g, tartaric acid 890 mg, citric acid anhydrous 720 mg, sodium citrate anhydrous 630 mg.

AUST 8 29693

8 sachets



Uses of citric acid:

Uses in the food industry:

> Used as acidulant in the manufacture of jellies, jams, sweet and soft drinks

> It is used for artificial flavoring in various foods

> Sodium citrate is employed in processed cheese manufacture

Uses of citric acid:

Uses in the cosmetic industry:

> It is used in astringent lotions such as aftershave lotion because of its low PH

> Citric acid is used in hair rinses and hair and wig setting fluids.

Procedure:

➢Weigh accurately (12.5g) of citric acid and dissolve it in (1L)of distilled water.

➤ and then take (3ml) of citric acid and titrate with
(0.25N) sodium hydroxide solution using 2 drops
of phenolphthalein solution as the indicator.

≻ Titrate until you get a pink color.



Calculation:

1. Calculate the quantity of citric acid present in our sample:

✓ When the base neutralizes the acid (at the endpoint), the number of equivalents of acid = the number of equivalents of base, therefore we can calculate the normality of citric acid from the following equation :-

```
(N1 * V1)Citric acid = (N2 * V2) NaOH
```

Where:

N1 :the normality of citric acid solution (unknown)

V1 :the volume of citric acid solution used (in our lab : 3ml)

N2 :the normality of NaOH solution (in our lab: 0.25)

V2 : volume of NaOH solution (the volume from coming out of the burette)

Calculation:

✓ Then, calculate the weight of citric acid from the following equation:

$$N = \frac{wt}{eq.wt} \times \frac{1000}{v}$$

Note:

Eq.wt of citric acid = 64.4g (why?)

Calculation:

2. Calculate the percentage (w/w) of our citric acid sample:

<u>note</u>: x sample weight

$$\% = \frac{wt \ of \ citric \ acid}{sample \ wt} * 100$$

Thank You For Listening