



Organic Pharmaceutical Chemistry

3rd Stage 2nd Semester

Lab No:1



The Preparation and Standardization of 0.1 N Potassium Permanganate Solution

Prepared by:

Assistant. lecturer: **Nadia Khalid**

Introduction:

- **Potassium permanganate** (KMnO_4 , m.wt.=158.0g/mol): is a dark purple or brownish black powder or almost black crystals. It is soluble in cold water and freely soluble in boiling water.
- Potassium permanganate is a **strong oxidizing agent** and It is decomposes on contact with certain organic substances.



Introduction:

- Due to its oxidizing abilities it has disinfectant and deodorizing properties. it is also astringent.
- It is widely used as a standard (volumetric) oxidizing solution because of its intense color which serves as an indicator in titrations besides its low cost.


Aim of the experiment:

- One liter of 0.1 N potassium permanganate solution is to be prepared.



Chemical principle:

- Standardization of potassium permanganate with anhydrous sodium oxalate (as the primary standard **why?**) follows **oxidation- reduction reaction**.
- In which potassium permanganate is the oxidizing agent where as sodium oxalate is the reducing agent.



Chemical principle:

- The titration is carried out in acid medium(why?).



Procedure:

A. Preparation of 1000 ml of 0.1 N KMnO₄ solution:

1- Weight accurately the appropriate amount of potassium permanganate (3.16g, **why?**).

$$N = \frac{Wt}{eq. wt} * \frac{1000}{V} \quad eq. wt = \frac{M. wt}{n}$$

$$0.1 = \frac{Wt}{\frac{158}{5}} * \frac{1000}{1000}$$



2- Transfer into a 250- mL beaker containing water and stir thoroughly breaking up the crystals with a glass rod, to effect solution.



-Add more water to the beaker and repeat the process several times until all the potassium permanganate has dissolved.



Procedure:

- Transfer the solution by using a funnel in to a one- liter volumetric flask, leaving the undissolved residues in the beaker.



5- Finally when the powdered KMnO_4 is complete dissolved, make up the solution up to the mark 1L(1000ml) by addition of water, then store in a dark clean closed- container



Procedure:

B. Standardization with Sodium Oxalate :

Weight accurately **0.1g** of anhydrous sodium oxalate (**M.wt=134g/mol**) previously dried to **110°C** and dissolve in **100ml** of water.

Transfer the solution and add sulphuric acid into a conical flask and then heat (about 70°C, **why?**) **Because the Oxidation of Sodium Oxalate is Rapid enough at such temperature.**


Procedure:

Add KMnO_4 from the burette with continuous stirring of the solution in the conical flask.

The first few drops result in a pink color persisting for about 20 seconds, wait until the color disappears and then continue the titration in the usual manner.

The end point is reached when a faint pink persists for about 30 seconds upon shaking the flask .

Record the volume of KMnO_4 solution used in the titration .



Calculation:

$$(N1 * V1) \text{ KMnO}_4 = (N2 * V2) \text{ Na}_2\text{C}_2\text{O}_4$$

Where:

N1 :the normality of KMnO_4 Solution to be Calculated (**unknown**)

V1 :the volume of KMnO_4 the volume descender from the burette

N2 :the normality of Sodium Oxalate solution (**in our lab 0.01 N**)

V2 : volume of Sodium Oxalate solution (**(in our lab : 5ml)**)

Thank
you