Surgical Instrument



Assis. Lect. Mohammed R. Abuljaleel

Surgical Instrument Classification

- •Cutting instruments: scalpel, scissors, rongeurs, and periosteal elevator
- •Bone instruments: chisels, osteotomes, gouges, saws, bonecutting forceps, curettes, and trephines
- •Grasping instruments : including needle holder, crushing and non crushing tissue forceps, hemostatic forceps, thumb forceps, towel clamps, and bone-holding forceps;
- •Retractors : hand-held and self-retaining
- •Suction tips
- Miscellaneous instruments

Cutting instruments

Scalpel contain two part

- **1. Handle** No.3 or No4
- 2. Blade
- ▶ No.3 = 10-15
- ▶ No.4= 20-25

Used for cutting tissue.



Figure 13-3 The four blades that fit a #3 scalpel handle (i.e., 10, 11, 12, and 15) and four blades that fit a #4 scalpel handle (i.e., 20, 21, 22, and 23). (Courtesy Miltex, Inc., York, PA.)



- •Scissors classified:
- Point type (sharp-sharp, sharp- blunt, and bluntblunt)
- Blade shape (straight and curved)
- Cutting edge type (plain and serrated)

Tissue or dissecting scissors



Mayo Scissors








General purpose scissors used for cutting tissue & material



Bandage scissors Used for cutting bandage



Suture scissor Used to removal suture material



Bone instrument

Bone catcher

Used to cutch the bone .





17. Small Reduction with Point Forceps



16.Small Reduction Forceps



14.Small Self-centering Forceps

Bone cutter




CR Used to cut bones



Bone drill Used to penetrate bones.

R Manual



R Electric




Bone plate and screw Used: fixed the fracture .

DCP Locking Plates, Cloverleaf Plate and Reconstruction Plate. Bone Screw, Cortical Screw, Cancellous Screw, Malleolar Screw DHS Screw, Interlocking Screw and Locking Screw.



Intra-medullary pin Used : to fixed fracture .



Retractors


Rand-held and self-retaining


***** Used to expose surgical site.

Hand-Held Retractors


Rake retractor Dever retractor


Different Retractor





Self Retractors Balfour retractor: used in abdominal surgery.





Finochietto retractor used : in thoracic surgery.



Grasping instruments

Needle holder
Simple N.H.
Compound N.H

- Spring N.H
- Used to catch needle during suturing

Needle Holder

Compound N.H.



Neck Inspect for wear, cracks, and missing inserts Inspect neck for cracks, blood, and stains Shanks


Mayo-Hegar Needle Holder

Inspect for wear, cracks, missing inserts Neck Inspect neck for cracks Test scissor Inspect for cracks,blood, and stains - Shanks Ratchets 11 Olsen Hegar Needleholder

Spring N.H




Tissue thumb forceps used to hold tissue



Dressing thumb forceps used to hold dressing and materials



Splinter thumb forceps used to remove the splint



Babcock grasping forceps used to hold bowel and bladder




Allis grasping forceps used to hold heavy fibrous tissue




Kocher grasping forceps To hold heavy, tough tissues




Thyroid tenculum (Lahey thy.for.) used to hold thyroid gland



Hemorrhoidal or lung clamp To hold hemorrhoid and lung



Sponge holding forceps To hold bandage or as retractor




Artery or hemostatic forceps Used to crush blood vessels



Doyen Intestinal forceps Used to hold intestine



To fixed the drapes



Gastric grasping forceps Used to grasp stomach



Suction tips

SUCTION


•Suction is used to remove blood and fluid from the surgical site. To reveal specific bleeding points for ligation or cautery.

•**Poole suction tips** have multiple openings and are best suited for use in body cavities.

• Frazier suction tips are often used during neurosurgical and selected soft tissue procedures (e.g., bulla osteotomy).

•The Yankauer suction tip may also be used






Miscellaneous instrument

Ruminal set
Ruminal magnetic
Ruminal magnetic
Trocar and cannula
Groove director
Groove director
Ear clamp
Kidney Shape
Surgical dram
Surgical trolley

Ruminal(weingart) set used in rumenotomy .



Ruminal magnetic

Used to catch the metallic material



Trocar and cannula Used: To emergency relief of free gas bloat






Use as guide to incision and suture





used to limit the line of incision of ear trimming





Surgical dram To save sterile instrument



Surgical trolley surgical instrument on it.





Thanks to Listen